```
1 // Dave Barry, Francis Crick Institute
 2 // 2018.01.17
  3 // david.barry@crick.ac.uk
 5 // Splits all files in input directory into constituent channels and saves results in separate sub-directories
          format = ".tif"; // specify image format
 9
         directory = getDirectory("Choose input files"); // get input directory
         print("Input: " + directory);
fileList = getFileList(directory); // get file list
print(fileList.length + " files.");
11
12
13
14
         outputDirectory = getDirectory("Choose output directory"); // get output directory
print("Output: " + outputDirectory);
15
17
18
          setBatchMode(true): // supresses windows opening
19
         for (i = 0; i < fileList.length; i++) {
    file = directory + File.separator + fileList[i];
    if(endsWith(file, format)){        // check if file is correct format
        print("\nFile " + file + " is a recognised format - processing.");
    run("Bio-Formats Importer", "open=[" + file + "] color_mode=Default rois_import=[ROI manager] view=Hyperstack"); // open image with bioformats
    getDimensions(width, height, sizeC, slices, frames); // get image dimensions
    print("Number of channels: " + sizeC);
}    if(ciseC = 1){</pre>
20
21
23
24
26
27
                     if(sizeC > 1){
28
                           run("Split Channels"); // split the image into constituent channels
29
30
                     subDirs = newArray(sizeC);
                     subDirs[c] = outputDirectory + "C_" + (c + 1);
if(!File.exists(subDirs[c])){
31
32
33
34
                                 File.makeDirectory(subDirs[c]);
35
36
37
                      titles = getList("image.titles"); // get the list of open images i.e the individual channels
                     for(j = 0; j < titles.length; j++){
    selectWindow(titles[j]);
    filename = subDirs[j] + File.separator + titles[j] + "_C" + (j + 1) + ".ome.tiff";</pre>
38
39
40
                           print("Saving " + filename);
run("Bio-Formats Exporter", "save=[" + filename + "] use compression=Uncompressed"); // save image
41
42
                           close();
43
44
               } else {
45
46
                     print("\nFile " + file + " is not a recognised format - skipping.");
47
48
          close("*");
49
50
51
          print("\nFinished");
          showStatus("Finished.");
setBatchMode(false);
52
53
```